





UroColor[™] Series

UroMeter 720[™]

UroMeter 120™





SD UroColor **Series

- Clear positive and negative results
- Good resistance to humidity
- Checking the interference in Glucose, Blood, Bilirubin, and Nitrite by Vitamin C ingested
- No color change after testing is completed
- No color spread between pH and Specific gravity
- No interference in protein by alkali urine
- Long shelf life: 24 months
- Different color packages by parameters
- Fast result visually or instrumentally



Product

Packing size : 50's, 100's, 150's, 500's, 1000's

• Shelf life : 24 months

UroColor Control

Positive and negative strip

Table of SD UroColor parameters

Cat. No.	Product	Blood	Bilirubin	Urobilinogen	Ketone	Protein	Nitrite	Glucose	рН	Specific gravity	Leukocyte	VitaminC
10UK01G	UroColor 1G							•				
10UK01P	UroColor 1P											
10UK01K	UroColor 1K				•							
10UK01B	UroColor 1Bi		•									
10UK01IU	UroColor 1IU			•								
10UK02	UroColor2					•		•				
10UK02K	UroColor 2K				•			•				
10UK03	UroColor3					•		•	•			
10UK03A	UroColor 3A	•				•		•				
10UK03K	UroColor 3K				•	•		•				
10UK03U	UroColor 3U			•		•		•				
10UK03N	UroColor 3N					•	•	•				
10UK04	UroColor4	•				•		•	•			
10UK04K	UroColor 4K	•			•	•		•				
10UK04S	UroColor 4S							•		•		
10UK04P	UroColor 4P				•	•		•	•			
10UK05	UroColor5	•		•		•		•				
10UK05B	UroColor 5Bi	•	•			•		•				
10UK05L	UroColor 5L	•				•	•	•			•	
10UK05K	UroColor 5K	•			•			•				
10UK05N	UroColor 5N	•				•	•	•	•			
10UK06	UroColor6	•			•	•	•	•	•			
10UK06B	UroColor 6Bi	•	•		•	•		•				
10UK06U	UroColor 6U	•		•	•			•				
10UK06L	UroColor 6L	•					•	•			•	
	UroColor 6NBi	•	•		•	•	•	•				
10UK07	UroColor7	•	•	•	•	•		•	•			
10UK08	UroColor8	•	•	•	•	•	•	•	•			
10UK08A	UroColor 8A	•	•	•	•	•	•	•				•
10UK09	UroColor9	•	•	•	•	•	•	•		•		
10UK09A	UroColor 9A	•	•	•	•	•	•	•	•			•
10UK10	UroColor 10	•	•	•	•	•	•	•	•	•	•	
10UK11	UroColor1 1	•	•	•	•	•	•	•	•	•	•	•

SD UroMeter 720™

- Auto calibration with power-on
- High capacity Max. 720 tests per hour
- Easy to input ID and use data
- Easy to input multiple ID by key board, PC and barcode reader
- Maximum 2,000 test results memory
- Easy software version up (PC ⇒ UroMeter 720)
 (No memory replacement)
- Applicable to your multi-language by MidiPC
- Easy to set-up & maintenance
- 4 ~ 11 parameter testing available

Fast result printing

- High speed thermal printer
- Complete data about patient and the result
- Highlightening abnormal result for quick review

LED indicator

- CIS module sensor
- Urine color & clarity determination : Eliminating the non-specific result from turbid urine

Various testing parameters

- Streamlining operations by immediate start-up
- It works with various urine chemistry dipsticks (parameters 4 ~ 11)



Specifications

Power Adapter : DC 12V 3A (80-250V, 50/60Hz)

• Size: 320 x 260 x 178mm (W x D x H)

Weight: 2.4 kg (5.3 lb)

■ Ambient operating temperature range (18-30°C)

Memory: 2,000 results

• Print: Internal line thermal printer

■ LCD: 8-line 21-character

RS232c Interface :

COM1: Communication with PC

COM2: Communication with barcode system



Various interface ports

- RS-232c, Parallel printer port
- Up/Downloading test result by Bidirectional interface
- A handheld bar-code reader identifying sample



- Bright and clear LCD
- 8-line by 21-character
- Automatic switching to testing mode from standby mode
- Easy to learn and operate



Continuous & automatic dipstick loading conveyor

- Take exact position through moving
- Easy maintenance by dismantling

SD UroMeter 120™

- Auto calibration with power-on
- Simple and excellent compatibility by USB communication port
- Enhancing user-convenience even for left-hander by ergonomics design
- Enhancing work efficiency by flexible options between quick mode and normal mode
- Easy to input multiple ID by key board, PC and barcode reader
- Applicable to your multi-Inguage by MidiPC
- 4 ~ 11 parameter testing available
- Maximum 300 tests/hour, Average 120 tests/hour
- Easy to learn and operate



- SD Urine Chemistry Analyzer
- Take increased throughput and improved workflow together with precision.

SD UroMeter 120™

Specifications

■ Power Adapter : DC 12V 3A (80-250V, 50/60Hz)

■ LCD: 8 line 21 character

■ Interface : RS232c :

COM1: Communication with PC

 ${\sf COM2}: {\sf Communication}$ with barcode system

USB Port





- 156-68 Hagal-dong, Giheung-ku, Yongin-si, Kyonggi-do, Korea 446-930
- TEL:+82-31-899-9700 FAX:+82-31-899-9740
- www.standardia.com sales@standardia.com